

CHIRANA T. Injecta, a.s.
Nám. Dr. Schweitzera 194
916 01 Stará Turá
Slovenská republika / Slovak Republic

IČO: **36794619**
IČ DPH/VAT No: **SK2022402822**

DECLARATION (no. EN 16-2020)

We, **CHIRANA T. Injecta a.s.**,

nám. Dr. A. Schweitzera 194
916 01 Stará Turá
Slovakia

manufacturer of:

STERILE SINGLE USE MEDICAL DEVICES

declare, that product


Sterile PERFUSION syringe - CHIRANA®

Types:

- 20 ml transparent for dialysis machine use (CH03020PT)
- 30 ml transparent for dialysis machine use (CH03030PT)
- 50 ml LL transparent (CH03050LL)
- 50 ml LL transparent with needle 2,0x30mm (CH030502030LL)
- 50 ml LL opaque (CH03050LLO)
- 50 ml LL opaque with needle 2,0x30mm (CH030502030LLO)
- 50 ml PERFUSION transparent (CH03050PT)
- 50 ml PERFUSION transparent with needle 2,0x30mm (CH030502030PT)
- 50 ml PERFUSION opaque (CH03050PO)
- 50 ml PERFUSION opaque with needle 2,0x30mm (CH030502030PO)
- 50 ml PERFUSION transparent BB (CH03050PTB)
- 50 ml PERFUSION transparent BB with needle 2,0x30mm (CH030502030PTB)
- 50 ml PERFUSION opaque BB (CH03050POB)
- 50 ml PERFUSION opaque BB with needle 2,0x30mm (CH030502030POB)
- 50 ml PERFUSION transparent F (CH03050PTF)
- 50 ml PERFUSION transparent F with needle 1,8x40mm (CH030501840PTF)
- 50 ml PERFUSION opaque F (CH03050POF)
- 50 ml PERFUSION opaque F with needle 1,8x40mm (CH030501840POF)

are suitable for use with infusion pumps defined in the table below.

in Stará Turá, date 24.04.2020


Rastislav Broska, PhD.
Quality Manager



Product type (REF code)	Setting of the infusion pump	Suitable infusion pumps	
20 ml transparent for dialysis machine use (CH03020PT)	20 ml BD Plastipak 20 ml Fresenius Heparin Syringe	ALARIS: GS, GH, GH plus, GH Guardrail, CC, CC plus, CC Guardrail, TIVA, PK, Aseña IVAC: 1000, 2000, 3000, 4000, 5000, 6000, 7000, P1000, P2000, P3000, PCA TERUMO: TE-311, TE-312, TE-331, TE-332, TE-371, TE-372 FRESENIUS: Pilot A2, Pilot C, Pilot Delta, Pilot Anaesthesia2, Orchestra Module DPS, Orchestra module DPS+, Injectomat Agilia, Injectomat MC Agilia, Injectomat TIVA Agilia B-Braun: PERFUSOR compact, PERFUSOR compact S, PERFUSOR Space MEDIMA: S1, S2, S300 FRESENIUS: 4008S (dialysis machine): 20 ml transparent for dialysis machine use (CH03020PT) 30 ml transparent for dialysis machine use (CH03030PT) ASCOR-MED: SEP 11S, SEP 21S, AP 12, AP 22, AP 14, AP 23, AP 24+	
30 ml transparent for dialysis machine use (CH03030PT)	30 ml Fresenius Medical Heparin		
50 ml LL transparent (CH03050LL)	50 ml CHIRANA original		
50 ml LL transparent with needle 2,0x30mm (CH030502030LL)			
50 ml LL opaque (CH03050LLO)			
50 ml LL opaque with needle 2,0x30mm (CH030502030LLO)			
50 ml PERFUSION transparent (CH03050PT)	50 ml B-Braun Omnifix		
50 ml PERFUSION transparent with needle 2,0x30mm (CH030502030PT)			
50 ml PERFUSION opaque (CH03050PO)			
50 ml PERFUSION opaque with needle 2,0x30mm (CH030502030PO)			
50 ml PERFUSION transparent BB (CH03050PTB)	50 ml B-Braun Original Perfusor Spritze	ALARIS: GS, GH, GH plus, GH Guardrail, CC, CC plus, CC Guardrail, TIVA, PK, Aseña IVAC: 1000, 2000, 3000, 4000, 5000, 6000, 7000, P1000, P2000, P3000, PCA TERUMO: TE-311, TE-312, TE-331, TE-332, TE-371, TE-372 FRESENIUS: Pilot A2, Pilot C, Pilot Delta, Pilot Anaesthesia2, Orchestra Module DPS, Orchestra module DPS+, Injectomat Agilia, Injectomat MC Agilia, Injectomat TIVA Agilia, Injectomat C-PS/CP-PS B-Braun: PERFUSOR compact, PERFUSOR compact S, PERFUSOR Space, PERFUSOR Secura FT, PERFUSOR F, PERFUSOR Secura FT µl, PERFUSOR Secura, PERFUSOR Secura P, PERFUSOR VII, PERFUSOR M, PERFUSOR Me, PERFUSOR M1, PERFUSOR ED, ED2, EDL, EDL2, PERFUSOR Bolustikolyse, PERFUSOR FM Clinico: Syritrol 12 MEDIMA: S1, S2, S300 ASCOR-MED: SEP 11S, SEP 21S, AP 12, AP 22, AP 14, AP 23, AP 24+	
50 ml PERFUSION transparent BB with needle 2,0x30mm (CH030502030PTB)			
50 ml PERFUSION opaque BB (CH03050POB)			
50 ml PERFUSION opaque BB with needle 2,0x30mm (CH030502030POB)			
50 ml PERFUSION transparent F (CH03050PTF)	50ml Fresenius Kabi AG Injectomat		ALARIS: GS, GH, GH plus, GH Guardrail, CC, CC plus, CC Guardrail, TIVA, PK, Aseña IVAC: 1000, 2000, 3000, 4000, 5000, 6000, 7000, P1000, P2000, P3000, PCA TERUMO: TE-311, TE-312, TE-331, TE-332, TE-371, TE-372 FRESENIUS: Pilot A2, Pilot C, Pilot Delta, Pilot Anaesthesia2, Orchestra Module DPS, Orchestra module DPS+, Injectomat Agilia, Injectomat MC Agilia, Injectomat TIVA Agilia, Injectomat S, Injectomat 50, Injectomat C-IS/CP-IS B-Braun: PERFUSOR compact, PERFUSOR compact S, PERFUSOR Space MEDIMA: S300
50 ml PERFUSION transparent F with needle 1,8x40mm (CH030501840PTF)			
50 ml PERFUSION opaque F (CH03050POF)			
50 ml PERFUSION opaque F with needle 1,8x40mm (CH030501840POF)			